

PREVENTIVE MAINTENANCE INSPECTION

racine

MEDICAL ELECTRONICS

148 Hill St. Shelton CT. 06484

Phone 800-836-0696

fax 203-924-7277

email racine.medical@snet.net

Name

Address

City

St.

Zip

Contact

SYSTEM TYPE	SERIAL NO.	SOFTWARE REV.

Probe Type	Serial No.	Peripherals	Serial No.

GENERAL INSPECTION

ACTION REQUIRED

<input type="checkbox"/> Verify proper operation of system, and user functions <input type="checkbox"/> Clean inside and outside of system <input type="checkbox"/> Verify casters / wheels roll freely, lock and are secure <input type="checkbox"/> Check operation of all fans <input type="checkbox"/> Clean camera optics and vcr heads, check operation <input type="checkbox"/> Perform all system diagnostics & gray scale test <input type="checkbox"/> Clean and check operation of B/W & color printers <input type="checkbox"/> Inspect internal and external cables for stress, discoloration, contact	<input style="width: 100%; height: 20px;" type="text"/> <input style="width: 100%; height: 20px;" type="text"/> <input style="width: 100%; height: 20px;" type="text"/> <input style="width: 100%; height: 20px;" type="text"/> <input style="width: 100%; height: 20px;" type="text"/> <input style="width: 100%; height: 20px;" type="text"/> <input style="width: 100%; height: 20px;" type="text"/>
--	---

ELECTRICAL INSPECTION

<input type="checkbox"/> Inspect AC power cord	<input type="checkbox"/> Measure AC leakage	<input type="checkbox"/> Verify AC on Aux. Equip.	<input type="checkbox"/> Measure AC
<input type="checkbox"/> Measure all DC voltages and ripple component			

IMAGE QUALITY

<input type="checkbox"/> Scan approved phantom, certify image uniformity, and system / transducer sensitivity	<input type="checkbox"/> Verify proper doppler & color operation
---	--

Comments

Field Engineer Customer

Date Date